


## ***Endocavity Color Doppler***

***USB&Wifi Handheld Untrasound***




**C10RQ**



<b>Frequency</b>	6.5-8MHz
<b>Depth</b>	30-120mm
<b>Mode</b>	B Mode B+M PW Doppler Color Doppler Power Doppler
<b>Applications</b>	Obstetrics / Gynecology



# Download our ultrasound software from the App Store or Google Play.



✓ Unlimited users

✓ Free software and unlimited software  
updates

✓ 18 months warranty

\*Extendable up to 3 years

## Software Controls

GN(gain)

D(depth)

ENH(enhancement)

DR(dynamic range)

F(frequency)

FocusPos

PRF

WF

Mode

8TGC

Biopsy

Annote

## Measurement

B: Length,Area,Circumference,GA  
(CRL,BPD,GS,FL,HC,AC)

B+M: Heart Rate, Time, Distance

B+PW:Velocity, Heart Rate(2), S/D, Depth



Package includes:

1\* wireless probe

1\* USB charger

1\* Type-C cable

1\* wireless charger

## Connection Type

*Type-C cable / WiFi (802.11n/20MHz/5G)*

*Probe built-in WiFi signal*

## Supporting Devices

*Tablet, Smart phone, PC, etc*

## Working system

*Android, IOS, Windows*

## Lithium Battery

*4200 mAh. 3.7V*

*Working time: ≥3h*

*Standby time: ≥12h*

*Battery life 5-8years*

*wireless charging*

## Dimensions

*345(h)x70(w)x30(d) mm*

Net weight      204g

Elements        192

Footprint        R10

## Boot-up

*Platform dependent, generally <5 seconds*



- \* **Training how to install the machine, training how to use the machine.**
- \* **Spare parts replacement and other software/hardware trouble shooting**
- \* **Offer OEM&ODM service for distributors**
- \* **Offer SDK depends on different models**
- \* **Free software and free upgrade the software**
- \* **18 Months Warranty, extendable up to 3 years**

**BEIJING KONTED MEDICAL TECHNOLOGY CO.,LTD**

**Tel:8610-60219113**

**Fax:8610-60219213**

**Service and Cooperation :info@konted.cn**

**Website:www.konted.cn**

**Facebook: Beijing Konted**

**Youtube: Medical Konted**

**Address: Room 111,1F, Building 3, No. 27 Yongwang Road,Daxing Biological Pharmaceutical Industry Base,Daxing District,Beijing 102629,China.**